

Inclusive Education

HOW-TO: Locking (and Unlocking) LAMP's Vocab and Menu Areas

LAMP Words for Life is customizable and is a huge benefit in providing a person with a communication system that is uniquely able to meet their vocabulary needs. There are many levels of security to lock in the edits (and general organization of the vocabulary), which protects the vocabulary from inadvertent changes. Here are instructions on how to navigate those locks.

*****It is important to consult with the speech-language pathologist about types of changes to be made to personalize the vocabulary.**

To Get Out of Guided Access (iPads' lock feature)

1. Triple click the home button
2. Enter your child's passcode (usually it is the school year, i.e. 202526).
3. Select END to get out of guided access

To Turn OFF an Editing Passcode in LAMP *(If your student usually has a passcode on editing on their talker, please remember to put it back on when you are done!)*

1. Choose Menu in the top right banner
2. Enter the passcode (likely the school year, i.e. 202526)
3. Choose Settings
4. Under Editing, turn Use Menu Passcode to OFF

To Turn ON an Editing Passcode in LAMP

1. Choose Menu in the top right banner
2. Choose Settings
3. Under Editing, turn Use Menu Passcode to ON

To Turn Guided Access (iPads' lock feature) Back On

1. While in LAMP, triple click the home button